

Restorative Materials

Expanded Functions Training
Lecture

Naval Dental Center, Pearl
Harbor

Materials we use...

- Amalgam
- composite resins
- cements
- porcelains
- RMGI, GI
- DBA's
- fluorides
- acrylics
- bleach
- Gypsum
- impression materials
- metals
- waxes

Amalgam

What's in it?

- Hg-mercury-dissolves solids, activates the reaction recombining the elements to make an alloy
- Ag- silver-reacts with Hg to form matrix; strength
- Sn-tin-low expansion but decreases strength
- Cu- copper- strength and hardness
- Zn-zinc- eliminates corrosion

Amalgam

- Admixed or Dispersed Phase
 - “sticks and balls”
 - mixture of lathe-cut and spherical
 - greater condensation force required
 - easiest to obtain contact in Class 2 preps
 - examples: Dispersalloy®, Valiant PhD®

Amalgam

- Spherical
 - “ball bearings”
 - easiest to condense, use largest condenser
 - “best” sealing margins
 - problematic in achieving tight contacts
 - examples: Tytin®, Valiant®

Can you mix different
amalgams in the same fill?

Yes!

- Large Build-ups are a good reason to mix amalgams
 - Start with fast set spherical for excellent marginal adaption and to allow timely removal of the matrix band, and
 - End with slow set admix to achieve tight contacts and allow time to carve

Amalgamation/ Trituration

Why is my amalgam too dry?

- Time: admix mixes more quickly than spherical
 - Under Triturated
 - Grainy & Dull
 - sticks to the capsule
 - Over Triturated
 - hot & crumbly
 - increases setting contraction, creep & decreases strength, working time

Pins

- Tooth crazing
- Regular TMS .031" diameter .027" drill
- Minim: .024" w/ .021" drill
- Minikin: .019" w/ .017" drill
- Minuta: .0145" w/ .0135" drill

Dentin Bonding Agents

Why do we use them?

Amalgams

- Seals margins
 - Amalgam shrinks .07%

Composites

- Micro and macro mechanical bonding

Dentin Bonding Agents

- Bifunctional Molecule
 - bonds to both “water loving” dentin and “water hating” composite
 - HEMA, 4 META
- Hybrid Layer: interpenetrating network of resin & collagen

Amalgam Bonding

- All Bond 2 ®
- Single Bond®
- One Step®
- Resinomer®

Composites

- Resin
 - Matrix
 - BIS-GMA: Highly viscous but very soft
 - absorbs water
 - Matrix diluent
 - TEGDMA: Decreased viscosity but increased filler load
 - Fillers--ground quartz
 - the more filler, the less matrix, the better the physical properties

Composites

- Coupling agents
- Polymerization Activator
- UV absorbers
- Inhibitors
- Pigments
 - Titanium and Aluminum Oxide

Composites

- Particle Size
 - Macrofill
 - Hybrid: e.g. Prodigy®, Herculite®
XRV
 - Microfill: e.g. Durafil®

Composites

- Clinical Factors
 - technique sensitive
 - isolation/moisture contamination
 - cavity prep-“C” factor
 - polymerization shrinkage
 - use with base or liner

Factors affecting the set of Composites

- Time
- intensity of light
- temperature
- distance
- resin thickness
- oxygen inhibition layer
- curing through tooth
- Shade of resin
- type of filler particle
- heat from curing light
- overhead fluorescent light

Microleakage

- Creates majority of the problems
 - post-op sensitivity
 - caries at margin
- Causes
 - hydrolysis of the bond
 - polymerization shrinkage
 - different (CTE: 30 vs 11) expansion rates

Glass Ionomers

- Glass--
 - powder
- Acid
 - --liquid

Glass Ionomers

- Conventional
- Resin Modified
 - light cured
 - Fuji II LC®
 - Vitrebond®

Sealants

- Moisture control
- cleaning tooth surface
- surface treatment
 - 37% etch 20 sec
 - 10 sec rinse & dry
 - bonding agent
- flow sealant into fissures
 - light cure min 20 sec
 - evaluate

Microabrasion

- 50 micron Al_2O_3 powder
- cleans preparation & fissures
- will create smear layer in dentin
- use for additional retention
- use to clean intaglio of crown

Bleach

- Hydrogen Peroxide: in office bleach
- Carbamide Peroxide: home bleach
- Sodium Perborate: walking tooth bleach
- Super Oxal 35% HP

Bleach

- Treat walking bleach cases with CaOH for 3 weeks before bonding composite

Fluorides

- 2% NaF
- 1.23%APF®
- Duraflor®
 - RD; Dry teeth; Apply (3X/ week); Let dry several minutes
- Gel-Kam® SnF
- Prevident® 1.1%NaF

Impression Trays

- Custom
- Stock
- Triple tray

Impression Materials

- Alginate: pour within 15 minutes
- VPS: pour after one hour < 7 days
 - Express®
- Polyether: pour < 7days; excellent moisture tolerance; stiffest material
 - Impregum ®
- Bite Registration
- Beware of bridges and bony exostoses

Hemostatic Liquids

- Hemodent--10% AlCl
 - Astringident--10% FeSO₄
 - Hemotrol
 - Visine--tetrahydrozoline
-
- use with knitted cords--2.5 X carrying capacity

Provisionals

- Jet® Acrylic
 - shrinkage
 - avoid moisture & contamination
 - strength
 - color stability