



U.S. AIR FORCE

Contagious Care Management



Maj Jeff Ellis
TSgt Jared Schultz
gt Benjamin L. Davis
16 OSS/OSM
30 Aug 04



U.S. AIR FORCE

Mission Capability Update

DRAFT FFGK8

Provides personnel and equipment for bare base operations in support of Special Operations Forces. Provides medical force protection insight of facilities at forward operating locations, particularly food, waste, water treatment, billeting, medical and sanitation facilities. Performs food and water source selection/evaluation, field hygiene/sanitation, vector/pest risk assessment, and initial occupational/environmental health site assessment for SOF personnel. Performs specialized small team/aircrew patient decontamination in far forward environments. Performs limited nuclear, biological, and chemical (NBC) hazard risk assessments. Adds Physician Extender support for clinical medicine activities as well as medical intelligence, medical administration, and operational planning, dependent on substituted AFSCs. The following substitutions are authorized: AFSC 4N051 with 041A3, 4A051, 4B051, 4E051, AFSC 042G3 with SEI 496 qualified 4N071, AFSC 4B071 with 043E3, and AFSC 4E071 with 043H3. UTC can deploy for 30 days without resupply. Deploys to MOP



U.S. AIR FORCE

Backup Slides



U.S. AIR FORCE











Treatment and Transport of CBRNE/TIC/TIM/DEW Exposed Personnel

Capability Description

Ref: AFSOC MAP 2033, Cap Shortfall #268

- Treatment (remote) of exposed personnel
- Transport personnel rapidly <1 min
 - Multiple (20) casualties at one time onto Fixed/Vertical takeoff aircraft
- Eliminate exposure to personnel and A/C

Operate Under any Env - OPS 2

- COA 1    
 - COA 2    
 - COA 3    
- (COA 1 + OTL Initiative)
- COA 3
(COA 2 + MF Initiatives)

Current Status/Risk

- OSM in process of purchasing 10 Man-Cages and 5 Casualty Care Systems
- 05 UFR submitted for \$50K for 10 CCS
- Current Risk - 4 (Catastrophic/seldom)
 - Personnel at Risk PAR - 20 patients
- Health Impact - 8 (severe, disabling) MFACs/Issues
- TTP Probability High
- Logistics. Who will maintain?
- Equipment disposal and decontamination
- COA if breach detected?



U.S. AIR FORCE

Multiple Casualty Isolation/Transport/Treatment

- Chemical Protective Barrier (Chem-Barrier) to contain chemical agent during transport of contaminated personnel and small equipment
- The personnel or small equipment are placed inside the Chem-Barrier. The Chem Barrier Door is closed and secured
- The Chem Barrier is designed to absorb liquid agents, thereby eliminating the contact hazard to personnel and equipment outside of the Chem Barrier. The vapors are contained in the pore structure of the Chem Barrier's fibers
- The Chem Barrier consists of an air permeable tent constructed of carbon cloth, butyl coated fabric and supports. The carbon cloth has been modified to enhance its durability, yet still retain the permeability needed for the personnel inside to breathe easily.





U.S. AIR FORCE


Single Casualty Isolation/Transport/Treatment



CASEVAC TRANSLOAD



You are here





CCS System Overview

U.S. AIR FORCE

- One-time use isolation unit to transport, hold, or protect casualties from contamination or cross-contamination.
- Provides protection from both chemical and biological warfare agents utilizing multiple reactive fabric laminates.
- Materials provide protection from contamination for 24 hours.
- Interior fabric is the made from same materials as outer fabric. Biological and chemical agents on the interior will also be absorbed for patient protection.
- One Medical Interface Bulkhead (MIB) containing two IV ports and one oxygen hose port.
- Filtered air is provided by two Micronel powered air purifying blower units using two standard NATO thread filter canisters and lithium batteries. Each blower provides four cubic feet per minute of air flow. Air flow can be altered to produce positive or negative air flow.



U.S. AIR FORCE

Photos - Side View

Right Side View





U.S. AIR FORCE

Photos - End Views

Head End



Foot End






U.S. AIR FORCE

Photos - Head End

Exterior - Head End



Interior - Head End





U.S. AIR FORCE

Photos - Head End

Exterior - Left Side of Head End



Exterior - Right Side of Head End







U.S. AIR FORCE

Photos - Foot End

Exterior - Foot End



Interior - Foot End





U.S. AIR FORCE

Photos - Foot End

Exterior - Left Side of Foot End



Exterior - Right Side of Foot End





U.S. AIR FORCE

Photos - Close-Ups

Medical Equipment Interfaces at
Head End



Medical Equipment Interfaces at
Foot End





U.S. AIR FORCE

Photos - Waist Section

Waist Section - Left Side





SOCOM Modifications

U.S. AIR FORCE

- Four sets of glove ports. One each located at head (for intubation), chest, waist, and feet.
- Light tan inner liner for faster identification of bleeders.
- Increased length and height for large framed patients.
- Woodland camouflage exterior.
- Disposable tear-off patient drape for protection of CCS and provider while loading contaminated patient.
- An extra generic equipment interface opening at head and foot. (Openings must be sealed with epoxy if used)



U.S. AIR FORCE

Recommended Changes

- Add accessory pockets to interior foot and head walls to store tubing, leads, and M-295/M-291 decontamination kits.
- Add Velcro straps to exterior foot and head walls to tie down tubing and leads.
- Create custom interface for ventilation tubing.
- Create custom interface for Propaq leads.
- Add another generic equipment interface for future communication system connections.



U.S. AIR FORCE

SOF FAWG Slides



Global Situational Awareness Tool (GSAT)

U.S. AIR FORCE

Capability Description

- Ref: AFSOC MAP 2033, Shortfall #275, 276
- Provide operational medical planners and civil engineers with a tool that can process and analyze existing databases and real-time data to assess ESOH risks.
- Assessment provides medical surveillance of war fighter exposures in-theater, relevant vaccines information, and prophylactic protocols that protect the health and safety of the war fighter

or follow-on technical assistance in hazard mitigation

	GOAS	Short	FYDP1	FYDP2
FYDP3				
COA 1 (Current COA)				
COA 2 (COA 1 + OTL Initiatives)				
COA 3 (COA 2 + MF Initiatives)				

Current Status/Risk



- SAIC contract - \$500K (disease layer, SIPR)
- SAIC contract - \$1M (Radiological Layer)
- Current Risk - 6 (Critical/Likely)
 - Personnel at Risk PAR - 90%
 - Health Impact - 8 (severe, disabling)
 - Probability - Seldom

LIMFACs/Issues

- Access to data and O&M cost to maintain
- Interface to other platforms (GEMS, DOEHRs, TOMS)



U.S. AIR FORCE

Tell GSAT what your missions are...




GSAT





**Then,
GSAT will tell you the hazards of
doing it**

[illegible]



Deployed Environmental Monitoring and Hazard Mitigation

U.S. AIR FORCE

Capability Description

- Ref: AFSOC MAP 2033, Cap Shortfall #284
- Robust, lean, agile field detection GPS locate, identification, and monitoring
 - TIC/TIM/CBRN in air/soil/water/food
- Ability to characterize hazards and provide user initial COA for hazard mitigation

Current Status/Risk

- Updating allowance std for FFGK8 to reflect new and improved capabilities
- 05 UFR submitted for \$320K
- Current Risk - 4 (Catastrophic/seldom)
 - Personnel at Risk PAR - 90%
 - Health Impact - 8 (severe, disabling)
- Probability - Seldom

Analysis Shortfall while DPAOR

- Ability to provide baseline for follow-on technical assistance in hazard mitigation
- COA 2 (COA 1+ OTL Initiatives)
- COA 3 (COA 2 + MF Initiatives)

LIMFACs/Issues

- Equipment integration requires multiple contractors working together



U.S. AIR FORCE

Env Monitoring & Hazard Mitigation

- Real-time device; Analyze liquids, solids, gases; user friendly, easy to maintain, Light portable, easily deployed into a far-forward environment



TODAY



Short Term 1-2 yrs



1-2 FYDPs



Sensing/
Detecting
Analyzing/
Reporting
Acting



HHA Kit



www.wasser1a.at



U.S. AIR FORCE

Process - Determine Capability Risk



		Probability				
		Frequent	Likely	Occasional	Seldom	Unlikely
Severity	Catastrophic	10	8	6	4	3
	Critical	8	6	4	3	2
	Moderate	6	4	3	2	1
	Negligible	3	3	2	1	1



U.S. AIR FORCE

Risk Definitions

(AFPAM 90-902, Operational Risk Management Tools)


- Probability of needing capability in future
 - Frequent: Continuously experienced (up to weekly)
 - Likely: Occurs regularly (weekly to annually)
 - Occasional: Occurs sporadically (2-5 times in 10 years)
 - Seldom: Rarely occurs (probable in 20 years)
 - Unlikely: Probably will not occur within 3 FYDP
- Severity if capability not improved
 - Health impact
 - Population at risk



Risk Definitions (Cont.)

(AFPAM 90-902, Operational Risk Management Tools)

U.S. AIR FORCE



Health Impact*

0 - None, annoyance, perceived quality
1-2 - Nuisance, irritant
3-5 - Temporary reversible illness or loss of capacity
6-7 - Permanent, non-severe illness or loss of capacity
8-10 - Permanent, severe, disabling, irreversible

* Adapted from AFI 91-301, Attachment 1

UNCLASSIFIED

























Integrity - Service - Excellence

DRAFT 26



U.S. AIR FORCE

Impact Based Upon Available Capability & Risk

<u>Capability</u>	<u>Current</u>	<u>FYDP</u>	<u>FYDP2</u>	<u>FYDP3</u>
Capability 1				
Capability 2				
Capability 3				
Capability 4				
Capability 5				
Capability N				

See Example